



## Bipolar resection electrode




Compatible to  
Olympusresectoscope

**TURIS**



item code	755-2
description	Loop electrode 90°
unit	12
Umax	1,8 kVp




WA22302D  
WA22503D

item code	755-1
description	Ball electrode
dimensions	Ø 3 mm
unit	12
Umax	1,8 kVp



WA22351C

Compatible to Olympus TURis  
Workingelements  
WA22366A + WA22367A



The resectoscope will be inserted through the urethra. With the TURis, diseased prostate tissue will be removed under sight. The resection succeeds with HF-loop electrodes which you plug onto the working element or the transporter in order to move it around. In a continuous steam the removed tissue layers will be directly transported out of the surgeon's sight by the assistance of a rinsing solution.

